

Aim :- To study the different properties of carbohydrates and determine an unknown carbohydrate.

Materials Required:

Chemicals

- | | |
|----------------------------|--|
| 1) Distilled Water | 7) Benedict's reagent |
| 2) Sodium Hydroxide (NaOH) | 8) Tollen's reagent |
| 3) Conc. Sulphuric acid | 9) Iodine reagent |
| 4) Iodine solution | 10) Water bath |
| 5) Conc. Nitric acid | 11) Test tubes and pipettes. |
| 6) Fehling's soln A & B | 12) Starch, Glucose, Sucrose
Lactose. |

Theory

- Carbohydrates are widely distributed in plants and animals. They have important structural and metabolic roles.

- Chemically carbohydrates are aldehyde or ketone derivatives of polyhydric alcohols.

- Glucose is the most important carbohydrate; the major metabolic fuel of mammals (except ruminants). And a universal fuel of fetus.

- It is the precursor for synthesis of all the other carbohydrates in the body.

Different types of Carbohydrates are :-

1) Monosaccharide - simplest form

2) Disaccharides - two monosaccharides

3) Oligosaccharides - three to ten monosaccharides.

4) Polysaccharides - condensation products of ^{many} monosaccharides


Teacher's Signature

Solubility Test



- Take a small amount of glucose, lactose, sucrose and starch in four test tubes A, B, C and D.
- Now add small volume of distilled water to the four test tubes and then shake the test tube well.
- We can see that glucose, lactose and sucrose are soluble in water whereas starch is insoluble in water.

Molisch's Test

- Take a small quantity of aqueous solutions of glucose, lactose, sucrose and suspension of starch in the four test tubes.
- Now add a few drops of Molisch's Reagent to the four test tubes.
- In the test tube A, pour conc. sulphuric acid slowly along side the test tube in small amount.
- In a similar way put conc. sulphuric acid of small amount in other test tubes as well.
- The Molisch's Reagent forms a coloured product that appears as a purple ring at the interface between the acid and the test layer.



Carbohydrate




5 - Hydroxymethyl
Furfural



Molisch Test



5 - (Hydroxymethyl)
Furfural



α - naphthol



Purple - Coloured
Dye.


Fehling's Test

- Take a small quantity of aqueous solutions of glucose, lactose, sucrose and suspension starch in the four test tubes A, B, C & D.
- Using a dropper add a small quantity of Fehling's solution A into the test tubes A, B, C and D.
- Now using a dropper, add a small quantity of Fehling's solution B into the test tubes A, B, C, D.
- Heat the test tube in a boiling water bath for sometimes.
- Glucose and Lactose form red precipitate of Cu_2O and others don't. Hence, Glucose & Lactose are reducing and others are non-reducing.


Benedict's Test

- Take a small quantity of aq. solution of the four carbohydrates in the respective test tubes.
- Add Benedict's Reagent to each Test Tube.
- Heat the test tubes in boiling water bath.
- Reducing sugars like Glucose and Lactose form red precipitate and others don't as they are non-reducing.

Fehling



Benedict's



Tollen's Test

- Take a small quantity of solution in each test tubes of the four sugars.
- Add a small amount of Tollen's Reagent.
- Heat the test tubes in boiling Water Bath.
- Glucose and Lactose produce silver mirror and no such silver mirror is produced in sucrose and starch.



Iodine Test

- Take a small quantity of aqueous solution of the four carbohydrates in respective test tubes.
- Add a small quantity of iodine in each test tube.
- Starch Reacts with iodine to form iodine-starch complex and others don't form any complex.

CONCLUSION




We come to know about a lot of properties of carbohydrates and test to differentiate reducing and non-reducing sugars.

Tollen's



Tollen

Reagent



Elemental
Silver

Answer